



# 3.1 Software Release Notes

For Navigator™ Software version 3.1

## Hardware Requirements

The software operates current WAVE 4500, WAVE 3500 and WAVE-MD instruments and runs on any computer that was running a previous version of the software.

## Software Requirements

The 3.1 client and server software can be installed on a new system or can upgrade version 2.1 or later.

The upgrades will not work automatically with versions earlier than 2.1 or with unpublished versions. Contact Transgenomic Customer Service for instructions on how to perform these upgrades.

**Note:** Navigator 3.1 will only be available for Windows XP machines. The last version of software that is available for Windows 2000 machines is Navigator 2.2. There will be no more updates for Windows 2000 machines.

## Upgrade Instructions

Two upgrade guides are now provided in PDF format. A *Navigator Server Upgrade Guide* should be included on the Server Upgrade CD. In addition, a new *Client Upgrade Guide* should be included on the client CD.

Please read all upgrade instructions carefully for 3.1

## New Features in Version 3.1

**Table 1: New Features**

Analysis	Added Sequence attribute to Analysis table.
Injections	Added Sequence attribute to Injection table.
Installation	The 4500 WAVE instrument can now be controlled on a Windows 7 computer.
Fragment Collector	Support for column and column S-pattern has been added.
Maintenance	Added ability to delete old parts.
Reporting	<ol style="list-style-type: none"> <li>1. Make Injection Report available from the analysis table.</li> <li>2. Added Resolution attribute to Analysis report.</li> </ol>

## Major Fixes in Version 3.1

**Table 2: Major Fixes**

Interface	<ol style="list-style-type: none"> <li>1. Stop injections if database files are full.</li> <li>2. Fix 3500 Fluorescence communication problems by implementing analog signal support.</li> </ol>
-----------	--

## Known Issues in Version 3.1

**Table 3: Known Issues**

DNA tab	<p>In some cases, imported web data (from MutationDiscovery.com) does not show features when the format button is pressed. As a workaround, please close the client and log back in.</p> <p>Deleting sequence across line breaks may inappropriately prevent the user from editing the end of the sequence. The user may select another sequence, then reselect the initial sequence to continue editing.</p>
---------	---

## Questions About 3.1

Please contact Transgenomic Technical Support by visiting: <http://www.trangenomic.com> and clicking on the Contact Us button on the menu bar.

Document Number: 482062-03